

Dental Instruments

Diagnostic

Instrument Image

Specification/Description




DE-01-000

Explorers

Double explorers comprised # 17 for gingival , calculus detection and #23 for carry detection

Use: For detecting and exploring abrasion



DE-01-001

Probes

Length of 6 mm for standard control point

Use : For checking the depth of periodontal pocket vertically or horizontally



DE-01-002

Explorer/Probes, 6mm

Standard control point for diagnostic patient



DE-01-003

Mirror Handle

Variety of mirror with size, style and handles, double mirrors, octagonal mirrors, mirror head front surface, anterior, posterior

Dental Instruments

Periodontal

Instrument Image

Specification/Description



Sickle/Jacquette Scalers

Anterior premolar/posterior, townner, white side, odont, morse

DE-01-004

Use: For removing calculus on the surface of gingival area



Black Titanium Instruments

Length of 6 mm for standard control point

DE-01-005



Sickle Scalers/Curettes

Facial/lingual, buccal/lingual

DE-01-006

Use : For removing supragingival calculus



Derbinement/Furcation

Buccal and lingual, mesial and distal, anterior, labial lingual

DE-01-007

Use :For removing residual calculus after ultrasonic scaling for deep pockets

Dental Instruments

Periodontal Surgery

Instrument Image

Specification/Description



DE-01-008

Surgical and Knives/files

Kirkland, Orban, crane Kaplan , sanders type– KM3/4, KCK-3, KOR1/2L

Use: For removal of soft tissue, and for removal of granulomatous



DE-01-009

Surgical Curettes

Kirkland, Kramer Nevins(KN3/4, KN4), sugarman(3/6,9/10, 11/12) Periotome

Use : For removal of remaining tissue of infected gingival area after extraction



DE-01-010

Periotomes

PT-1, PT-2, PT-3, PT-4 (titanium regular, titanium angle)

Use : For removal of teeth with minimum tissue damage



DE-01-011

Surgical Chisels

13k/TG– 3.5/3mm(3.5/4.5mm), 13k/TG titanium, KN1/2, C36/37, SOLT 3/4

Use :For scraping out and reshaping bone; for removal of flaps and splitting bone

Dental Instruments

Implant instruments

Instrument Image



DE-01-012

Specification/Description

Osteotomes

Concave STR, Convex CVD
(2mm,3mm,3.3mm,4mm,5mm),

Use: For bone spreading for implant



DE-01-013

Sinus Lift instruments

Bonescoop 5mm, periosteal ,bone plugger
3x1.5mm, bone graft 1.0/1.7mm (titanium,
serrated)

Use : For separation of sinus lining from bone



DE-01-014

Osteotomes 2

OSTM(3mm-4.7mm), OSTM STR
(1.8mm-4.7mm), straight, Jovanovic

Use : used in surgical procedures for preparation
of implants



DE-01-015

Endo Forceps

Endo Forceps (75mm, 95mm, 110mm)

Use :For implantation

Dental Instruments

Implant instruments

Instrument Image

Specification/Description



Bone Syringe

Bone syringe (3mm,3.5mm,4mm)

Use: For implantation

DE-01-016



Castroviejo

Castroviejo STR (8.5cm,17cm)

Use : For implantaion

DE-01-017



Tissue Plier

Adson (12cm, 15cm), Crane Kaplane 15cm,
self locking 16cm, tissue forceps T/C

Gerald 18cm

Use : For preventing slipping while reducing the
weight over 20%

DE-01-018



Tissue, Cotton and Dressing pliers

Atraumatic 18cm, Allison baby 12.5cm,
Allison 15cm

Use :For implantation

DE-01-019

Dental Instruments

Implant instruments

Instrument Image



DE-01-020

Specification/Description

Cotton Pliers

Adson (12cm, 15cm), adson baby of 12.5size

Use : For implantation support



DE-01-021

Titanium instruments

Titanium Molt-9, PR-3, 24G, Buser

Use : For implantaion



DE-01-022

Accu Laxating

Accu lux (straight, curved, spade, distal, mesial, large distal, large mesial)

Use: For pulling , lifting up the periosteum and loosening the tooth from the lip



DE-01-023

Scalpel

Scalpel of standard 180 degree and 16.2cm

Use : Forr implantation

Dental Instruments

Operative

Instrument Image

Specification/Description



PKT Waxing

Gold Fogel –1,3R

Use: For shaping, carving or pouring molten wax

DE-01-024



Amalgam Carriers

Small 1.5mm, medium 2mm, large 2.5mm, Jumbo 3mm

Use : For carrying and dispensing amalgam materials

DE-01-025



Gingival retractors/cord packers

GRM-(2,3,4), GR-10, cord packers serrated

Use : For protecting tissue during cavity preparation and placement of composite fillings

DE-01-026



Plastic filling

Anterior /Posterior—darby perry trimmer (3R, 21L,49, H2, B2, IPCOA)

Use :For composite placing and shaping

DE-01-027

Dental Instruments

Endodontic

Instrument Image



DE-01-028

Specification/Description

Pluggers and condensers

RCP set of (0.3, 0.4, 0.5, 0.55, 0.6, 0.65, 0.75)
for anterior / posterior–
8P, 8.5P, 9P, 9.5P, 10P, 10.5P, 11P, 11.5P

Use: For vertical pressing and compacting filling materials



DE-01-029

Spreaders and heat carriers

RCS sets(RCS 30, 40, 50, 60), heat carriers
(0.20, 0.25, 0.35, 0.55)

Use : For filling solid medicine materials into the endodontic area during lateral condensation



DE-01-030

Rubber Dam Clamps

Anterior (6, 9, 9s, 212),
Premolar (0, 1, 1A, 2, W2, 2A, 2AS, 27N, 29), lower molar (3, W3, 7, W7, 7A), upper molar (56, 56S, 4,

Use : stabilizes and holds the dam secure



DE-01-031

Ivory Light weight

Ivory light weight 17cm, brewer 17.5cm, stoke 17.5,

Dental Instruments

Surgical

Instrument Image

Specification/Description



DE-01-032

Periosteals

Bennet, P-149 , titanium molt-9, PR-3, P15,P16,24G, HF-20, periostal elevater

Use: For retracting periosteum after removing gingival tissue



DE-01-033

Surgical Curettes (Spoon)

86 lucas, 87 lucas, 87 lucas, (8-12)miller, molt-2,molt-4,molt -5L,molt-5R

Use : For curettage and cyst removal



DE-01-034

Bone Chisels/files

Single bevel BF1x(B1 bevel,miller), Gardner (1,2,5,5R,6,9), Chandler (1,2,4) , (235g,130g) mallet, mead-2mallet

Use : For final smoothening/removing bone by pushing/pulling out



DE-01-035

Chisels and Gouges

Partsch 17cm, freer 16cm, lucas 16cm

Use :For scooping bone

Dental Instruments

Surgical

Instrument Image

Specification/Description



DE-01-036

Bone Rongeurs

Blumental 15cm, blumental 45degree, 14.5cm , blumental 90degree, 14.5cm

Use: For tooth extraction as a whole by trimming and recontouring alveolar bone



DE-01-037

Root Splinter

Root fragments straight (45, 90 degree), peet sharp (15, 45, 90 degree)

Use : For extracting root tips of the tooth where normal forceps cannot reach



DE-01-038

Sponge and Dressing Forceps

With out rathet (14.5cm,18cm,20cm), curved (17.5cm, 21cm)



DE-01-039

Haemostatic forceps

Pean 12.5, baby-crile(pean) 14cm, crile 14cm, Halsted 12.5cm, micro-halsted 12.5cm (curved)

Use :For lockable forceps for seizing the end of a blood vessel to control hemorrhage

Dental Instruments

Surgical

Instrument Image



DE-01-040

Specification/Description

Needle holder

Needle holders crile-murry 15cm, 18cm, Adson 17.5cm, Castraoviejo 16cm, Boynton 12cm

Use: For holding and guiding the needle securely during suturing



DE-01-041

Gum Scissors

Dean 16.5cm, Matzenbaum 14cm, locklin 16cm, mayo 14.5cm CVD, Lagrange CVD 11.5cm, suture 14cm, Castroviejo 10cm

Use : For cutting threads and gum



DE-01-042

Tungsten Carbide insert

Metzenbaum 14.5cm, goldman 13cm, wire cutting 12cm, dean 16.5cm , locklin 16.5cm, Kelly 16cm

Use : For cutting threads and gum



DE-01-043

Raspatories

Koning 15.5cm , williger 14cm , 16cm; tessier 20cm, 22cm,

Use : For hooking skin during surgery, also as retractors for holding back the lips and cheeks

Dental Instruments

Surgical

Instrument Image



DE-01-044

Specification/Description

Skin hook and retractors

Frazier 13cm, Joseph 16.5cm, Gillies 18cm, Standell-Stillle 18cm, Langenbech-Green 16cm, Roux 16.5cm

Use: For hooking skin during surgery, for holding back the lips, cheeks and tongue



DE-01-045

Mouth Gag

Mouth Gag(adult size) 14cm, child size 12cm, pedo size 10cm

Use : For keeping the mouth open



DE-01-046

Scalpel handles

Scalpel of size 13.5cm, 12cm, 14.5cm, double head 14.5cm, scalpel 180degree;16.2cm



DE-01-047

Use : For holding a mess blade

Syringes

Syringe 1.8ml, 2.2ml; paroject

Use : For injecting anesthetics from a cartridge; for supplying air, water or a mist in three ways

Dental Instruments

Surgical

Instrument Image



DE-01-048

Specification/Description

Elevators

Crane pick 8, seldin 1L, seldin 1R, seldin 4L, cryer M25, apical 3L, flohr sharp end

Use: For pulling, lifting up the periosteum and to loosen the tooth from the periosteal



DE-01-049

Root tip picks/teasers

Root 78 str, root 79L, root 80R



DE-01-050

Use : For removing residual root tips of tooth

Extraction Forceps

Lower molars cow-horn beak, lower molar left side, upper anterior 1AF, Upper molars hook handle 10H, 15, 16S, upper molars/incisors, root canines and premolars strong beak

Use : For extracting teeth from alveolar bone



DE-01-051

Extraction Forceps/pedodontic pattern

158 upper molars universal, 160 lower molars universal, lower central and roots, upper central, 5 lower incisors

Use : For extracting teeth from alveolar bone

Dental Instruments

Crown instruments

Instrument Image



DE-01-052

Specification/Description

Crown Scissors

Straight saw edge 10.5cm, curved saw edge 10.5cm , universal 12cm

Use: For removing crowns permanently by cutting off the seal between tooth and crown after splitting



DE-01-053

Crown Spreader

Crown spreader 14cm, flat; teeth

Use : For splitting and spreading the crown open for removal



DE-01-054

Crown remover

Crown remover of size between 15cm and 20cm is of standard



DE-01-055

Use : For gripping onto the crown and removing from the tooth

Impression trays

Solid partial set of 10, perforated partial set of 10 (upper/lower), solid pull denture , perforated full denture

Use : For forming a crown frame

Dental Instruments

Orthodontic

Instrument Image

Specification/Description



DE-01-056

Cutters

Universal distal end cutter, distal end cutter slim/handle, wire cutter(hard,soft)

Use: For cutting wire during orthodontic procedures



DE-01-057

Band remover

Ticker/straight/long/angled/adhesive/bracket remover

Use : For removing orthodontic bands from teeth



DE-01-058

Pliers
Universal 16cm, Aderer 14.5, Double serrated 14cm, arrow clasp forming 12cm, pointed jaw 14cm, mosquito plier 12.5cm



DE-01-059

Use : For bending wire during orthodontic procedures

Matrix retainers

Ivory 8(5mm, 6mm,7mm), siqveland wide, tofflemire universal

Use : For holding metal bands or strips around the tooth creating a temporary mold.

Instrument Image



DE-01-060



Wax knives

Sizes of 12.5cm 17cm

DE-01-061



Use : For carving wax in laboratories

DE-01-062



DE-01-063